

Computational Fluid Dynamics simulations and in vitro testing on the importance of expiratory valves for spacers

Emeline Furon¹, Faustine Vanel¹, Thierry Porée¹, Jean-Christophe Dubus²

¹Laboratoire Protec'Som; ²Service de pneumologie pédiatrique, Unité de Médecine Infantile, CHU de Marseille - Hôpital de la Timone

Introduction

A valved holding chamber (VHC) is more than a spacer extending the distance between the mouth of the patient and the pressurized metered dose inhaler (pMDI). Fitted with one-way inspiratory valves, preventing exhalation into the device, some also present an expiratory valve¹.

Drug delivery from holding chambers is dependent on multiple factors, including the efficacy of the valve system². Using computational fluid dynamics (CFD) modelling and *in vitro* experimentation, the influence of 2 small apertures (3*3mm) situated before the inspiratory valve, as a mean of expiratory valve, was tested.

CFD modelling

A simple 3D geometry (figure 1) of a pMDI coupled to a VHC with a volume of 200 mL and comprising 2 opposite apertures (3 mm diameter) was created with SOLIDWORKS (DassaultSystèmes).




Figure 1 : Model VHC 3D geometry.

An automatic mesh generator (Protec'Som) for spacers was used to create a tridimensional mesh of the spacer. CFD modelling using a laminar flow was performed with the open source OpenFOAM software (OpenCFD Ltd) and the SimpleFOAM Solver.

CFD modelling predicted that, with a constant inspiratory flow of 15L/min, 69.3% of the air going through the mouthpiece would come from the side apertures, leaving only 4.6 L/min of the flow coming from the back of the spacer.

Conclusion

The *in vitro* validation of the computational modelling simulation of flow through a model VHC, confirms the pertinence of our *in silico* models.

Effective seal between the patient and inspiratory valve is essential, therefore an efficient valve system is required, especially for children with lower inspiratory flow.

Computational Fluid Dynamics modelling




Figure 2 : Side view of the model spacer and MDI: vertical velocity.

VHC Seal

Chamber seal and aperture influence was determined by placing the model spacer (with or without apertures plugged) in line between 2 flow meters and connecting to a pump (figure 3) delivering a constant flow (15 or 50mL/min). CFD modelling data were confirmed *in vitro*, where 66.8% of the airflow came from the apertures so only 33.2% coming from the back of the spacer. That proportion of air traversing the VHC increased to 53% when using a 50L/min flow.




Figure 3 : VHC seal testing experimental set up.




Figure 4 : VHC seal was determined by measuring the flow ratio at 15 and 50 L/min.

In vitro deposition

Lung deposition of fluticasone (Flixotide 50µg/dose, GSK) was assessed *in vitro* by connecting the model spacer to a model trachea and collecting filter (figure 5) coupled to a breathing simulator (Copley scientific) generating a child breathing profile (155mL/5 breaths). The presence of those apertures reduced the lung dose deposition *in vitro* by 38% (3.54 ± 0.27 vs 5.71 ± 1.46 µg/puff, $p < 0.005$) when compared to the spacer with plugged apertures.




Figure 5 : *In vitro* filter deposition experimental set up.




Figure 6 : Lung and throat deposition *in vitro* estimation

References

- Nikander, K., Nicholls, C., Denyer, J. & Pritchard, J. The evolution of spacers and valved holding chambers. *J Aerosol Med Pulm Drug Deliv* 27 Suppl 1, S4-23 (2014).
- Mitchell, J. P. & Nagel, M. W. Valved holding chambers (VHCs) for use with pressurised metered-dose inhalers (pMDIs): a review of causes of inconsistent medication delivery. *Primary Care Respiratory Journal* 16, 207–214 (2007).



21. ISAM-Congress
Santa Fe, NM
03rd - 07th June 2017

