

combiHaler®



*Spacer for mechanical
ventilation circuit*

combiHaler®

Combi-Haler® is a spacer designed to be used in invasive or non invasive mechanical ventilation circuits.

Combi-Haler® allows for the use of both a vibrating mesh nebulizer and a pMDI simultaneously.

Intended for adults suffering from Chronic Obstructive Pulmonary Disease (COPD), asthma or chronic obstructive respiratory failure, Combi-Haler® can be used **in hospital and for home hospitalization.**



Performances



Security



Simplicity



Practicity



Performances

Pulmonary deposition x 2 with a nebulizer and a pMDI (in comparison with a traditional assembly).



Simplicity

Combi-Haler® can be used for many hours.

With Combi-Haler®'s design and airtight caps, it can remain on the ventilation circuit for several hours.



Security

Transparency : Combi-Haler®'s transparent material allows for the monitoring of the circuit.

No change of the ventilator settings : reduces the risk of errors.

Single patient use : limits the risk of contaminations.



Practicity

Combi-Haler® was designed with a light material to not put weight on the ventilation circuit.

Thanks to its specific adaptor, Combi-Haler® can be used with a soft mist MDI (Respimat®).



French quality production

- ISO 13485 : certification of every processus of manufacturing and packaging
- 100% french products

combiHaler®

Use with nebulizer



(compared to Aerogen® T-adapter)



Combi-Haler® **double** the pulmonary deposition of Amikacin (with a nebulizer).

Use with pMDI

(compared to a T-piece)



Combi-Haler® **double** the pulmonary deposition of medication (with a pMDI).

Combi-Haler® reduces deposition by impaction and/or by sedimentation of substances in the spacer and the circuit.

Protec'Som SAS

Espace d'Armanville
24 rue du Train Renard - BP 46
50700 Valognes - France
Tél. : +33 (0) 233 887 094
Fax. : +33 (0) 233 887 095



Protec'Som Belgium

c/o CFCIB, avenue des Arts, 8
1210 Brussels - Belgium
Tél. : +32 (0) 78 050 017

Protec'Som North America Inc.

1416 Jean-Berchmans-Michaud St.
Drummondville, Québec, Canada J2C 7V3
Tél. : +1 819 850-2202
Fax. : +1 819 477-5512



@ProtecSom
<https://twitter.com/ProtecSom>



@LaboProtecSom
<https://facebook.com/LaboProtecSom>



www.laboprotecsom.com

